
GRAINMAN

EQUIPMENT No.64 LABORATORY SAMPLE SHELLER FOR PADDY RICE

*Model 64-115-60-WDC
Laboratory Paddy Rice Sheller
Shown with Dust Collector*

Designed to comply with service requirements of The U.S.D.A. Federal/State Rice Inspection Service, the Grainman No. 64 Laboratory Sample Sheller for Paddy Rice provides a rapid, easy and efficient method of shelling, while utilizing only one piece of equipment that will accommodate dry paddy rice samples.

It can be carried anywhere, requiring small space and may be plugged into any convenient electrical outlet.

The Grainman No. 64 Laboratory Sample Sheller for Paddy Rice is offered in two basic models: with dust collector and without dust collector. The dust collector recovers hulls in a clear plastic cup connected to the machine by two flexible hoses. The hulls material can be weighted later for other analysis. Also, the collector is recommended for cleanliness of operation but may be by-passed to blow the hulls away, if desired.

Our Equipment No. 64 greatly improves its well-known use to:

- Grain Elevators.
- Rice Mills.
- Grain & Seed Control Official Inspections.
- Universities
- Rice Seed Producers
- Rice Research Centers

All information above is subject to change without previous notice

GRAINMAN


OPERATION:

The rice paddy sample is placed in the hopper. From the hopper the paddy is feed continuously by the vibrating feeder between the rolls.

The fast steel roll has corrugations and the slow synthetic rubber roll has a smooth surface. The shear effect created by the fast and slow rolls is applied to the paddy kernel in the gap between the rolls thus tearing the husk apart. Kernels and husks are then separated by a controlled air force

A cyclone collector (-WDC Models) separates the husk from the air. The husks are deposited in a clear plastic container. If the collector is not used or available (-WODC Models), then husks can be blown away by mean of the fan.

FLOW DIAGRAM:

IMPORTANT FEATURES:

- **Shelling Efficiency from 95% To 100%:** Our No. 64 Sheller can give the highest shelling efficiency among all known methods used in the market. This is achieved by careful manufacturing and superb craftsmanship.
- **Only One Control for Any Paddy Variety:** By turning one single control, the gap between the two rolls is adjusted to accept any paddy variety, whether long grain (thin kernel) or short grain (thick kernel).
- **Shell any Dry Paddy Sample for Easy and Fast Inspection:** Make fast sample evaluation to speed up grading, Other methods of shelling paddy samples are impractical, Instant detection of objectionable kernels.
- **Seed Producers - Seed Users:** Our No. 64 Sheller improve your crop yield and grade by rapidly inspecting the seed to be planted. It provides the fastest and best method to determine the quality of a lot of seed rice.

MODEL No.	DESCRIPTION	OPERATION
64-115-60-WDC	Grainman No. 64 Laboratory Paddy Sheller <u>with</u> Dust Collector	115 Volts 60 Hz 1-phase
64-220-50-WDC		220 Volts 50 Hz 1-phase
64-115-60-WODC	Grainman No. 64 Laboratory Paddy Sheller <u>without</u> Dust Collector	115 Volts 60 Hz 1-phase.
64-220-50-WODC		220 Volts 50 Hz 1-phase

Shipping Information:

Part No.	Description	Weight	Dimensions (in)
64BOX-W	Wooden Box	200 Lbs	34" x 20" x 32"

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