GRAINMAN

EQUIPMENT No.61 RICE SIZING DEVICE

Model 61-115-60 Rice Sizing Device shown with two tray plates (tray plates sold separately)

Designed to comply with service requirements of The U.S.D.A. Federal/State Rice Inspection Service, the Grainman No. 61 Rice Sizing Device provides a rapid, easy and efficient method of broken separation from whole rice grains.

It requires a small space and may be plugged into any convenient electrical outlet.

The GRAINMAN® Rice Sizing Device accurately tests samples in order to determine the percentage of broken and whole rice. The *recommended sample* size is 1000 grams of cleaned shelled, milled rice.

This machine consist of a sturdy wood frame that accommodates perforated aluminum plates mounted on four metal *adjustable* pedestals and connected to an eccentric linkage, activated by an electric motor to produce a sharp *reciprocating lateral movement*. All firmly secured to a wooden base for table-top operation. Two plates are used in the machine simultaneously to make a sample run (sold separately).

Our Equipment No. 61 greatly improves its well-known use to:

- Grain Elevators.
- Rice Mills.
- Grain & Seed Control Official Inspections.
- Universities
- Rice Seed Producers
- Rice Research Centers



All information above is subject to change without previous notice

GRAINMAN

OPERATION:

1. Insert two perforated plates into the machine with the perforation size according to the type of separation required.
2. Weight 1000 grams sample of cleaned, shelled and milled rice.
3. Start the motor (be sure the plates are securely fastened before starting).
4. Place the rice sample into the top plate and place a sample pan under the hopper to catch the rice that flows over the plates.
5. The machine is run until the sample stops flowing over the plates into the sample pan. Stop the motor.
6. Remove and empty the rice that remains on the plates into a separate container. Change to the next set of plates and repeat the above steps, until all necessary separations have been made.
7. Take all the obtained portions (from plates and sample pan) and weight them.
8. Because the total sample was 1000 grams, the percentage of each portion is obtained by weighing it and simply moving the decimal point one place. The percentage of whole rice is taken from the weight of the rice which remains after the last set of plates have been used. From these figures you will obtain the percentage of separation.

MAIN PARTS:

Plates mounted on wood frame (sold separately)

NOTES:

- Plates are perforated with round holes, centers staggered. Wood frame mounted. Available Plates Sizes: 5/64", 6/64", 10/64", 11/64", 12/64", 0.135", 13/64",
- It comes with a TEFC electric motor with 3-prong H.D. chord.

| MODEL No. | DESCRIPTION | OPERATION |
|-----------|--|-------------------------|
| 61-115-60 | Grainman No. 61 Laboratory 1000 grams Rice Sizing Device | 115 Volts 60 Hz 1-phase |
| 61-220-50 | Grainman No. 61 Laboratory 1000 grams Rice Sizing Device | 220 Volts 50 Hz 1-phase |

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